

# AIA-600II Assay Specifications

## ST CPEP II Test Code 040

---

Screen	Item	Data
<b>Screen 1</b>		
	Lot	***Enter current cal. lot no.
	Cal 1	0 ng/mL
	Cal 2	0.5 ng/mL (example)
	Cal 3	2.0 ng/mL (example)
	Cal 4	6.0 ng/mL (example)
	Cal 5	15.0 ng/mL (example)
	Cal 6	33.0 ng/mL (example)
 <b>Screen 2</b>		
	Name	#CPR2
	Unit	ng/mL
	Smpl	10
	Dil	110
	2Reag	0
	Code	0
	Assay Range Low	0.04
	Assay Range High	30.00
	Reference Range Low	0.00
	Reference Range High	0.00
	DP (No. of decimal points)	2
 <b>Screen 3</b>		
	Code	40
	No. (Calibrators)	6
	Mul. (Replicates)	3
	Equ	13
	CV (Calibration curve stability)	90
	STAT (Analyte status)	0
	PRCL (Assay Protocol)	1
	Factor 1	
	A	1.000
	B	0.000
	Factor 2	
	A	1.000
	B	0.000
 <b>Screen 4</b>		
	<b>Dilution Factors:</b>	
	SP1 (Specimen 1)	1
	SP2 (Specimen 2)	10
	CAL	1
	CTRL	1
	CODE (SDS code)	40
	PR (Dilution mode)	3