

Tosoh Bioscience Inc

AIA-2000 Automated Immunoassay Analyzer

The AIA-2000 Automated Immunoassay Analyzer is the latest advancement in high throughput immunoassay testing from Tosoh Bioscience, South San Francisco, Calif. Tosoh has added new functionality to the AIA-2000 to increase workflow and minimize time. The user can load up to 960 tests (48 trays) in a new, easy to load hybrid sorter; increasing walk-away time to approximately 4 hours. Tip capacity has been increased to 576, and can be added or

changed at any time at the user's request. Sample racks can be loaded with any combination of 200 bar-coded primary tubes or sample cups.

The AIA-2000 features fully automated daily maintenance and record keeping, automatic inventory management and automated pre-treatment (for assays that require this step). In

addition, the analyzer features a 17-inch touch screen panel and improved user interface, allowing easy access to most of the controls from the screen. Multiple robotic arms, new staging platforms, and mechanical adjustments speed the overall sample processing.

As with all AIA systems, Tosoh's reagent test cup technology produces accurate and consistent results, requiring no on-board



the new marketplace for products, services, and career opportunities



**Electronic Imaging
Materials, Inc.**

BARCODE LABELS

**Adhere to all types of slides and tubes.
Our labels can also withstand extreme
chemical exposure.**

- * Blank labels for most printers
- * Pre-printed labels for immediate use
- * Single labels & multi-part sets
- * Complete printing system packages including, printers, scanners & software



For more information,
or to request a
label sample, please log
on to our website at
www.eiminc.com/lab
or call us at
800-535-6987



**Providing
Laboratory
Labeling
Solutions
Since
1987**

refrigeration, no pre-mixing, and no pre-measuring.

AIA-2000 system features include: throughput of 200 tests/hour, 576 tips (6 tip racks), 2 substrate bottles, walk away time of 4 hours, 3 incubators with different processing steps, 2 sampling arms for faster processing of second reagent, and a full test menu.

Tosoh Bioscience Inc

(800) 248-6764; www.tosohbioscience.us