

High-throughput automated immunoassay analyzer

The *AIA-2000 Automated Immunoassay Analyzer* allows users to load up to 960 tests (48 trays) in an easy-to-load hybrid sorter — increasing walkaway time to approximately four hours. Features include throughput of 200 tests/hour, 576 tips (six tip racks), two substrate bottles, three incubators with different processing steps, two sampling arms for faster processing of a second reagent, and a full test menu. The system includes fully automated daily maintenance and record keeping, automatic inventory management and automated pre-treatment for assays that require this step. The analyzer features a 17-inch touch-screen panel to allow for easy access to most of the controls.



Tosoh Bioscience

Visit www.rsleads.com/910ml-171