

Symmetry in Diagnostics



A New Day for Vitamin D



Fortify your Lab with Vitamin D Testing

Tosoh Bioscience is proud to introduce the ST AIA-PACK 25-OH Vitamin D assay for use on Tosoh AIA System Analyzers. Vitamin D measurement is intended as an aid in the determination of Vitamin D sufficiency.

The antibodies have equimolar affinity for Vitamin D₂ and D₃ and the assay measures Total Vitamin D. Dry reagent format ensures 90-day calibration stability for minimal waste and cost effective testing.

TOSOH BIOSCIENCE

www.tosohbioscience.us



Fortify your Lab with Vitamin D Testing

- The Tosoh ST AIA-PACK 25-OH Vitamin D correlates well with gold standard methods and measures 25-OH Vitamin D₂ and D₃ in equimolar proportions.
- 90 day calibration stability makes it feasible to run this test even with low volumes.
- Available on a broad range of analyzers – AIA-600II, AIA-900, AIA-1800 and AIA-2000.
- Utilizes Tosoh’s Unit Dose Test Cup technology.
- Single, unitized test cups require no pre-mixing, no pre-measuring and no on-board refrigeration.
- Test cups are bar-coded for easy identification and inventory management.
- Assay range of 4 to 120 ng/mL.
- Assay time of 40 - 50 minutes.
- Sample type - Serum, Na heparinized plasma or EDTA plasma.
- Low CVs.



AIA-900



AIA-900
with 9 Tray Sorter



AIA-900
with 19 Tray Sorter



AIA-2000

Part #	Product Name	Description
025234	ST AIA-PACK 25-OH Vitamin D	20 cups x 5 trays (100 tests)
025334	ST AIA-PACK 25-OH Vitamin D CALIBRATOR SET	Six levels, 2 x 6 x 1 mL
025734	ST AIA-PACK 25-OH Vitamin D PRETREATMENT SET	PR-1-1 x 32 mL, 6 Reagent vials (empty) PR-2-6 x 5 mL
025434	AIA-PACK 25-OH Vitamin D CONTROL SET	Two Levels, 2 x 2 x 1 mL



Tosoh Bioscience, Inc.

6000 Shoreline Court, Suite 101
 South San Francisco, CA 94080
 Tel: (800) 248-6764
 Fax: (650) 615-0415
www.tosohbioscience.us