

# Comprehensive Laboratory Solution for Hemoglobin A1c Testing



The 501RP+ Reporting Software, delivers workflow efficiency for your laboratory from processing to reporting. The HLC®-723G8, an automated glycohemoglobin analyzer, provides high precision quality performance.

Together, the instrument's reliability and accuracy with Tosoh's latest comprehensive data management system provide caregivers the confidence necessary to deliver the accurate results.

## PRECISION & ACCURACY

Paperless, high resolution chromatograms graphs with clear separation and less than 2% CVs

## SPEED & EFFICIENCY

Direct determination of SA1c in 1.6 minutes through ion-exchange, HPLC methodology

## TRACEABILITY

Barcoded positive patient identification and full reagent traceability and on-board stability



HLC®-723G8



501RP+ Reporting Software

For more information, call 1.800.248.6764

# TOSOH BIOSCIENCE

[www.tosohbioscience.us](http://www.tosohbioscience.us)

"Tosoh Bioscience" is a registered trademarks of Tosoh Corporation  
Part # ADS-HLC-001 Rev. 1/0619

