



## Accurate SHBG Results with Test Cup Convenience

*Quality Results  
for Quality Life*



Using total testosterone and SHBG, results can be provided as free androgen index (FAI), free testosterone\*, weakly bound testosterone\* and bioavailable testosterone\*.

*\*Requires albumin measurement*

**TOSOH BIOSCIENCE**

# Available for Labs and Workloads of All Sizes



**AIA-360**  
36 Tests Per Hour



**AIA-900 Benchtop**  
90 Tests Per Hour



**AIA-900 Loader Model**  
90 Tests Per Hour




**AIA-900 with 9 tray sorter**  
90 Tests Per Hour  
Also available with 19 tray sorter



**AIA-2000**  
200 Tests Per Hour

## Unit Dose Test Cup Reagents

- Single, unitized cups use a dry reagent format that requires no pre-mixing, no pre-measuring, no on-board refrigeration and no waste
- Dry reagent format ensures calibration stability of up to 90 days for most assays
- Reagents have a one year shelf life from date of manufacture
- The AIA test cups and trays are labeled with code and lot number for automated scheduling and inventory management on the AIA-2000 and AIA-900 sorter models
- ST AIA-PACK reagents provide a fast turnaround time due to a 10 minute incubation protocol
- Extensive test menu provides a wide range of immunoassay tests including thyroid, cardiac markers, tumor markers, reproductive hormones, anemia markers, metabolic markers, diabetes, kidney markers and more



**TOSOH BIOSCIENCE**

Tosoh Bioscience, Inc.  
3600 Gantz Rd  
Grove City, OH 43123  
Tel: (800) 248-6764  
Fax: (650) 615-0415

[www.tosohbioscience.us](http://www.tosohbioscience.us)