





Quality Immunoassays for Quality of Life

Clinical Testing Solutions with Consistency, Stability and Harmony

Unit Dose Test Cup Reagent System

- Dry reagent format
- 90-day calibration stability for most assays
- Single, unitized cups - No pre-mixing or measuring
- Interchangeable reagents between Tosoh AIA® Systems
- Bar-coded for identification/inventory management
- 10 minute incubation time for most assays
- Biotin free immunoassay designs - No risk of biotin interference
- Extensive test menu includes assays for cardiac, thyroid, reproductive, tumor, metabolic, diabetes, kidney, anemia, Vitamin D and more



TOSOH BIOSCIENCE

Tosoh AIA Test Menu

10-minute incubation unless otherwise noted.

Tumor Markers

AFP
CEA
PA
27.29
CA125
CA19-9

Diabetes

C-Peptide II
IRI

Thyroid

TSH
TT3
T4
FT3
FT4
T-U
TSH 3rd-Gen ♦
TPOAb ♦
TgAb ♦

Cardiac Markers

BNP
CK-MB
cTnl 2nd-Gen
Myoglobin

Anemia

FER
B12 ♦
FOLATE ♦
RBC FOLATE ♦

Reproductive

βHCG
HCG
DHEA-S
E2
FSH
LH II
PRL
PROG II
PROG III
SHBG
Testosterone

Kidney Markers

BMG
Cystatin C
intact PTH

Metabolic

CORT
HGH

Additional

ACTH
Homocysteine
IgE II
PAP
25-OH Vitamin D ♦

♦ > 10 minute time for pretreatment and/or incubation. Not available on the AIA-360.

One reagent for all Tosoh AIA Systems

The AIA-PACK® Test Cup format works with every Tosoh AIA Automated Immunoassay System. Laboratories with more than one Tosoh system can use the same reagents interchangeably between instruments. Transition between systems is seamless, ensuring consistent results and efficient, economical operation for the laboratory.



AIA-360



AIA-900
Benchtop



AIA-900
Loader Model



AIA-900
with 9 Tray Sorter
(Also available in 19 tray
sorter option)



AIA-2000

Not all immunoassays are the same. Using the right immunoassay gives you confidence in the results you deliver.

Choose wisely. Choose Tosoh.



TOSOH BIOSCIENCE

Tosoh Bioscience, Inc.
6000 Shoreline Court, Suite 101
South San Francisco, CA 94080
Tel: (800) 248-6764
Fax: (650) 615-0415
www.tosohbioscience.us

"Tosoh Bioscience", "AIA-PACK" and "AIA" are the registered trademarks of Tosoh Corporation